


PRECICE[®] HPRT Assay Kit Hypoxanthine-guanine phosphoribosyl transferase Assay Kit

PRECICE[®] HPRT Assay kit provides a simple non-radioactive assay for measurement of Hypoxanthine-guanine phosphoribosyltransferase (HPRT) activity in a convenient 96-well plate format. In the assay, HPRT activity is measured as a rate of production of IMP, which is oxidized by recombinant IMP-dehydrogenase with concomitant formation of NADH₂. The formation of NADH₂ is continuously monitored spectrophotometrically at 340nm. The assay is linear up to 340nmol/hour/ml of HPRT activity, enables the detection of as low as 6.75nmol/hour/ml and provides results with uncertainties of 6.5%. Specific activities of HPRT in lysates of erythrocytes and cultured cells measured by PRECICE[®] HPRT Assay kit are comparable to those obtained previously by radiochemical procedure.


Principle

IMP-dehydrogenase enzyme catalyzes the irreversible oxidation of IMP to XMP with simultaneous reduction of NAD to NADH₂. **PRECICE[®] HPRT Assay Kit** takes advantage of this property for monitoring the activity of HPRT. In the assay, HPRT activity is measured as a rate of production of IMP, which is immediately oxidized to XMP by IMPDH with concomitant reduction of NAD to NADH₂ measurable by absorbance at 340nm.



Accurate estimation of HPRT activity with in **NOVOCIB's** IMPDH-based assay is done:

- in unidirectional reaction where IMP produced is immediately converted by IMPDH to XMP;
- through continuous monitoring of HPRT activity at 340nm in a convenient 96-well plate;
- by simultaneous analysis of 15 samples in triplicate per plate (2h);



Calibration curve of IMPDH-based PRECICE[®] HPRT Assay

The rate of the increase in absorbance at 340nm per hour as a function of the concentration of human recombinant HPRT enzyme (NovoCIB, ref. E-Nov-9). The changes in absorbance corresponding to 3 different control hemolysates diluted in complete reaction mixture to final hemoglobin concentration 1mg/ml (n=4) are shown by filled squares, filled triangles and filled circles. The insert shows a linear correlation observed in whole range of 21ng/ml up to 1.5µg/ml of recombinant HPRT; the units are the same.

- Kit Content:**
- Bacterial IMPDH
 - IMPDH Reconstitution Buffer
 - NAD 100x
 - DTT
 - HKM Buffer containing hypoxanthine (10x)